


Tibia & Femur ^^ Proximal Tibia High Osteotomy Anatomic Plate

Proximal Tibia High Osteotomy Anatomic Plates are indicated for osteotomies, treatment of bone and joint deformities, fixation of fractures, and malalignment caused by injury or disease, such as osteoarthritis, of the


Generally seen after high-energy traumas in young people and low-energy traumas in osteoporotic bone in the elderly.

Proximal Tibia High Osteotomy Anatomic Plate, it is designed to rearrange the stress distribution by changing the direction of the load axis passing through the diseased knee joint area, which is subjected to continuous loading, and to stabilize the fractures and deformities. Two options wedge-free Osteotomy Wedges: 6mm – 8mm – 10mm - 12mm- 14mm




Proximal Tibia High Osteotomy Anatomic Plate Screws Info

Anatomically contoured, eliminating the need for intraoperative contouring and minimizing soft tissue irritation.



Elongated Combi holes in the neck and shaft facilitate plate adjustment and allow locking or compression options.



Tapered ends allowing submuscular plate insertion.



Proximal Tibia High Osteotomy Anatomic Plate Screws

Osteotomy Wedge

6 mm

8 mm

10 mm

12 mm

14 mm

5.0 mm Non-Locking Cortical Screw



5.0 mm Locking Cannulated Cortical Screw



5.1 mm Locking Cannulated Cancellous Screw Cortical Screw



5.1 mm Non-Locking Cancellous Screw

